
Name:	Date:
-------	-------

Continuous Renal Replacement Therapy



Across

- **4.** monitors for air bubbles within return line
- **7.** used to prevent the blood from clotting
- **8.** two bean shaped organs
- **9.** continuous renal replacement therapy
- 10. waste removal bag containing ultrafiltrate

- **11.** the functional unit of the kidney
- 12. the movement of solutes through a semi-permeable membrane from higher concentration to lower concentration

<u>Down</u>

1. monitors effluent line for presence of red blood cells

- 2. pressure sensor that measures pressure required to push blood back to the patient
- **3.** most common type of anticoagulant
- **5.** sudden loss of kidney function
- **6.** Dialysis fluid that is administered on the outside on the filter and directly enters the effluent bag