

Name: _____

Date: _____

Lab

M A C R O M O R P H O L O G Y B M B U G Q Y F I
I F J M O A T X P Y V A V T S C Q R W O P D I N
W V N V W P I V D G R K Q S V R A I U L S X I E
E L F T L X V D G N I X I F T A E H A H Y A P S
A L S W I S S K E K H C D G Y Z G N B O T T N W
O A Y L V L H U T M C G T B X N C C B S T E T E
P Z L Q G S T R S N L I D E J G L V M X Y L K N
D O D L C H Z R L W A A U V M E V A H U F O Q R
R Z O P Y B L C A L G T I Q Z A R N T M P I U I
G I H O I Z T G N O K L S T F G L G V E X V K C
Z Q T S G X K E T W V B N I N F B F Q S M L P H
S K E Z I F X Z S G G C Q A S E I X Z K P A H E
X U M L P Y V H C X P I N U E E R D Q C N T Y D
Z I T W I B C T U Z O F S K Q V R E C J Q S X M
J G N A L T D W N R Q W Z F M Z I Z F G M Y Q E
B L A D E C O L O R I Z A T I O N B U F T R Q D
D F R Y B Q T P J E I F O F R Z C X B U I C Z I
P Z D L J F Z A E S U E F Y B F Y H B A C D Q A
W P A E K Y D C B Q W Z E F K S O E C X P W K A
Y A U A Q R D C F G M N E D A P F A S B X N R B
U U Q F W U R Y Z Q E T K C V N V J K C Z C V U
N R L C F M E U F O A P E L B I T P E C S U S M
G J A B Z E Q C A L M D X V Z O Z O Q C Q W Z L
N C X R P G O H P K J D Z O I R P C P R A L S N

differential media
Decolorization
susceptible
Diff Quick
plate

quadrant method
Crystal violet
heat fixing
resistant
flame

macromorphology
enriched media
gram stain
Slants
Tube